

Characterizing Anatomical Variability And Alzheimer's Disease Related Cortical Thinning in the Medial Temporal Lobe

Long Xie, Laura Wisse, Sandhitsu Das, Ranjit Ittyerah, Jiancong Wang, David Wolk, Paul Yushkevich


PENN Image Computing and Science Lab (PICSL), University of Pennsylvania

ShapeMI Workshop - MICCAI 2018

Granada, Spain, Sept 20th 2018





Subregions of the medial temporal lobe (MTL)



- BA35 and BA36 are subregions of the perirhinal cortex (PRC)




Early neurofibrillary tangle pathology (NFT) and related cell/synapse loss of Alzheimer's disease begin in the medial temporal lobe (MTL)

Stages of Neurofibrillary Tangle Pathology [Braak & Braak 1991,95]







- Structural measurement of the MTL substructures are promising biomarkers of AD.

Hard to quantify due to the existence of discrete anatomical variants, defined by the folding and branching patterns of the CS



Angevine et al., The Human Brain in Photographs and Diagrams

Borders and extents of BA35 and BA36 depend on the depth of CS



- Failing to account for the anatomical variability in the analysis degrades our ability to reliably localize and accurately quantify brain regions in individual subjects.

Outline/Aims


- Automatic segmentation pipeline to segment ERC, BA35/36 from structural MRI
 - Apply to both T1-weighted and T2-weighted MRI
- The multi-template thickness analysis pipeline
 - To extract regional thickness of these structures
 - Establish anatomical meaningful correspondence between subjects
- Characterizing anatomical variability and Alzheimer's disease related cortical thinning
 - Apply the proposed pipeline to a large dataset of the baseline T1-weighted MRI scans from Alzheimer's Disease Neuroimaging Initiative (ADNI)

Outline/Aims


- **Automatic segmentation pipeline to segment ERC, BA35/36 from structural MRI**
 - Apply to both T1-weighted and T2-weighted MRI
- The multi-template thickness analysis pipeline
 - To extract regional thickness of these structures
 - Establish anatomical meaningful correspondence between subjects
- Characterizing anatomical variability and Alzheimer's disease related cortical thinning
 - Apply the proposed pipeline to a large dataset of the baseline T1-weighted MRI scans from Alzheimer's Disease Neuroimaging Initiative (ADNI)

ASHS: A multi-atlas segmentation pipeline optimized for MTL substructure segmentation in T2-weighted MRI

3 Tesla T1-weighted MRI ($1.0 \times 1.0 \times 1.0\text{mm}^3$)



3 Tesla T2-weighted MRI ($0.4 \times 0.4 \times 2.0\text{mm}^3$)



➤ Segmenting a new subject

● Inputs:

- High-resolution T2w MRI
- 1 mm^3 isotropic T1w MRI

● Algorithms:

- ANTS deformable registration²
- Joint label fusion³
- Corrective learning⁴


● Output:

- Segmentation of hippocampal subfields and MTL cortex in the space of the T2w MRI


➤ Atlas set: 29 subjects (15 controls, 14 MCI)

- High-resolution T2w MRI
- 1 mm^3 isotropic T1w MRI
- Manual segmentation in the T2w MRI space

ASHS can reliably segment MTL substructures





Generalized DSC: 0.779




ASHS-T1: From T2-weighted to T1-weighted MRI

- **Motivation: T1w MRI is the most commonly acquired MRI modality**
 - Although T1w MRI does not provide enough contrast to visualize hippocampal subfields, the MTL cortex can be reliably segmented in T1w MRI.
 - More than 1000 T1w MRI scans of subjects at different stages of AD are available in the Alzheimer's Disease Neuroimaging Initiative (ADNI) database
 - Large sample size allows us to the characterize cortical thinning patterns of anatomical variant
- **Difficulty: Dura has similar intensity as gray matter in T1-weighted MRI**




The dura is segmented as gray matter by FreeSurfer¹




ASHS-T1¹: A multi-atlas segmentation pipeline optimized for MTL substructure segmentation in T1-weighted MRI

3 Tesla T1-weighted MRI (1.0 x 1.0 x 1.0 mm³)




Upsampled T1-weighted MRI (0.5x0.5x1 mm³) and manual segmentations



- **Atlas set: 29 subjects (15 controls, 14 MCI)**
 - 1 mm³ isotropic T1w MRI
 - Manual segmentation in the upsampled 0.5x0.5x1 mm³ T1w MRI space
- **Segmenting a new subject**
 - **Inputs:**
 - 1 mm³ isotropic T1w MRI
 - **Algorithms:**
 - ANTS deformable registration²
 - Joint label fusion³
 - Corrective learning⁴
 - **Output:**
 - Segmentation of the MTL cortex in the space of the upsampled T1w MRI


Segmentation accuracy is comparable to that in T2w MRI



Dice coefficients between automatic and manual segmentations in SR-T1w and T2w MRI

Modality	ERC	BA35	BA36	Dura
T1w MRI	0.76 (0.03)	0.70 (0.06)	0.78 (0.04)	0.75 (0.05)


Dura mater is often labeled as gray matter by FreeSurfer¹



% of dura voxels in manual segmentation labeled as			
Method	Dura	Gray Matter	Other
FreeSurfer ¹	N/A	62.4 (10.5)	37.6 (10.5)

Extended to segment the anterior/posterior hippocampus and parahippocampal gyrus and it is publicly available

Generalized Dice Similarity Coefficient = 0.8357




<https://sites.google.com/view/ashs-dox/>

Outline/Aims


- Automatic segmentation pipeline to segment ERC, BA35/36 from structural MRI
 - Apply to both T1-weighted and T2-weighted MRI
- **The multi-template thickness analysis pipeline**
 - **To extract regional thickness of these structures**
 - **Establish anatomical meaningful correspondence between subjects**
- Characterizing anatomical variability and Alzheimer's disease related cortical thinning
 - Apply the proposed pipeline to a large dataset of the baseline T1-weighted MRI scans from Alzheimer's Disease Neuroimaging Initiative (ADNI)

Build template for each variant from the atlas set using graph-based groupwise registration¹




- Manual segmentations from the right side were flipped, yielding 58 samples
- Subtype of each manual segmentation was manually assigned

Build template for each variant from the atlas set using graph-based groupwise registration¹




- Pairwise highly regularized coarse registrations were performed within each subtype
- Edge weight was set to the GDSC of BA35, BA36 and CS after registration

Build template for each variant from the atlas set using graph-based groupwise registration¹



- Minimum spanning tree² was build for each subtype
- The sample that is closest to all the other samples was identified as the “root”

Drawback: the template is very similar to the “root”



“Root”

Template derived from
graph-based registration


Shape correction and model shape variability using pointset geodesic shooting

- Nice properties of large deformable diffeomorphic metric mapping (LDDMM) via geodesic shooting^{1,2,3}
 - The deformation field between template and the target subject can be compactly represented by the initial momentum
 - The space of initial momentum provides a linear representation of the high-dimensional non-linear diffeomorphic transformation space, which is important for statistical analysis of shape variability




<https://www.youtube.com/watch?v=jurwTuCJn7s>


Shape correction via pointset geodesic shooting¹




Shape correction via pointset geodesic shooting¹




Shape correction via pointset geodesic shooting¹



Shape correction via pointset geodesic shooting¹



The corrected template does seem to represent the mean shape




“Root”

**Template derived from
graph-based registration**


Shape corrected template

Principal component analysis on the initial momenta¹




Principal modes capture expected variability – T1 Atlas

Variant with continuous CS




Variant with discontinuous CS



- Shape analysis using the manual segmentations of the T1 atlas

Fit the template to a new automatic segmentation along the minimum spanning tree

- Automatic segmentation using ASHS-T1
- Register to all the manual segmentations in the atlas set
- Pick the most similar 6 atlases and perform a weighted vote to decide the variant membership of the new sample
- Connect the sample with the most similar manual segmentation in its group
- Warp the template to the target sample




Outline/Aims

- Automatic segmentation pipeline to segment ERC, BA35/36 from structural MRI
 - Apply to both T1-weighted and T2-weighted MRI
- The multi-template thickness analysis pipeline
 - To extract regional thickness of these structures
 - Establish anatomical meaningful correspondence between subjects
- **Characterizing anatomical variability and Alzheimer's disease related cortical thinning**
 - Apply the proposed pipeline to a large dataset of the baseline T1-weighted MRI scans from Alzheimer's Disease Neuroimaging Initiative (ADNI)

Applied to baseline T1-weighted MRI scans from ADNI

	Control	Preclinical AD	Early Prodromal AD	Late Prodromal AD	Dementia
N	180	94	139	109	81
Age (yrs)	72.0 (6.0)	74.5 (5.7) ***	73.0 (6.9)	71.7 (6.8)	74.9 (7.8) ***
Gender (M/F)	94 / 86	32 / 62 **	80 / 59	57 / 52	47 / 34
Education (yrs)	16.9 (2.4)	16.1 (2.7) *	15.7 (2.9) ***	16.6 (2.6)	15.4 (2.6) **
MMSE	29.0 (1.3)	29.0 (1.1)	28.0 (1.7) ***	27.2 (1.9) ***	23.2 (2.1) ***



AD-related cortical thinning – ADNI cross-sectional dataset



t-statistical maps of the contrast between amyloid-negative cognitively normal adults and the other 4 groups (preclinical AD, early and late Prodromal AD, dementia) of both variant-templates, with age and education as covariates. Only results of the right hemisphere are shown. The patterns are similar on the left side.


Early neurofibrillary tangle pathology (NFT) and related cell/synapse loss begin in the medial temporal lobe (MTL)

Stages of Neurofibrillary Tangle Pathology [Braak & Braak 1991,95]




- Being able to replicate findings from pathology studies supports the utility of structural measurement of MTL cortex in tracking early Alzheimer's disease progression

Classification using the initial momenta




Effect of AD on MTL shape – ADNI Dataset

Variant with continuous CS



Variant with discontinuous CS



- AD is associated with decrease in overall size of the MTL, cortical thinning and widening of the CS
- These shape features may provide complementary information in identifying disease groups

Summary

- The first automatic segmentation pipeline of MTL cortex using T1-weight MRI that explicitly accounts for the confound of the dura mater¹
- A novel multi-template analysis pipeline to quantify shape variability of anatomical variants of the MTL
- Progression of cortical thinning that is consistent with known progression of NFT pathology within the MTL cortex related to AD
- Proposed method may have important utility in the early detection and monitoring of AD and the findings in this study may help us better understand the effect of AD on the shape of MTL substructures

THANK YOU!

- **Link to ASHS and ASHS-T1 software**

<https://sites.google.com/view/ashs-dox/>

- **Publicly available atlas sets for T1w and T2w MRI**

<https://www.nitrc.org/projects/ashs/>

- **Acknowledgements**


NIH grant numbers R01-AG056014, R01-AG040271, P30-AG010124, R01-EB017255, R01-AG055005 and the donors of Alzheimer's Disease Research, a program of the BrightFocus Foundation



SCAN US!


T1 Atlas: derived from manual segmentation in T2 space

T1w MRI (1 mm³)




Super-
resolution
Upsample


SR T1w MRI (0.5x0.5x1 mm³)



Affine Alignment




Super-
resolution
Upsample



T2w MRI (0.4x0.4x2 mm³)

SR T2w MRI (0.4x0.4x1 mm³)




T1 Atlas: derived from manual segmentation in T2 space


T1w MRI (1 mm³)




SR T1w MRI (0.5x0.5x1 mm³)




Propagate labels



T2w MRI (0.4x0.4x2 mm³)



SR T2w MRI (0.4x0.4x1 mm³)





T1 Atlas: derived from manual segmentation in T2 space


T1w MRI (1 mm³)



SR T1w MRI (0.5x0.5x1 mm³)



T2w MRI (0.4x0.4x2 mm³)



SR T2w MRI (0.4x0.4x1 mm³)

